


INNOVATION
MADE BY GEUDER®

NEW

PURE EFFICIENCY PHACO

Greatest Efficiency for
Smallest Incisions

SET FOR COAXIAL MINI PHACO



More gentle insertion due to the smooth transition of the sleeve to the tip

Increased followability and efficient holdability due to improved fluidics

Short ultrasound times with minimal stress to the capsular bag due to easier entry into hard lens nuclei by the smaller edge geometry

Maximized active surface for ultrasound emulsification due to a 40° angle and a 3-step design

Geuder®
Precision made in Germany

PURE EFFICIENCY PHACO

Greatest Efficiency for Smallest Incisions

for 2.2 mm incision



G-24310 PEP 2.2 MM ULTRASONIC TIP
Pure Efficiency Phaco
for 2.2 mm incisions, titanium, 40°, 3 step design, purple



G-24311 PEP 2.2 MM INFUSION SLEEVE
Pure Efficiency Phaco
for 2.2 mm incisions, silicone



G-34088 SINGLE-USE 2.2 MM PHACO KNIFE
Pure Efficiency Phaco
angled, 2.2 mm, double bevel
6 pcs. per box, sterile



G-34176 SINGLE-USE 19G MVR STILETTO
double bevel, 19 Gauge / 1.1 mm
6 pcs. per box, sterile



G-34181 SINGLE-USE 20G MVR STILETTO
double bevel, 20 Gauge / 0.9 mm
6 pcs. per box, sterile



G-24318 PEP 2.2 MM I/A HANDPIECE
Pure Efficiency Phaco
for 2.2 mm incisions, straight design
with G-24311 PEP infusion sleeve



G-24418 PEP 2.2 MM I/A HANDPIECE
Pure Efficiency Phaco
for 2.2 mm incisions, angled design
with G-24311 PEP infusion sleeve

for 2.4 mm incision



G-24313 PEP 2.4 MM ULTRASONIC TIP
Pure Efficiency Phaco
for 2.4 mm incisions, titanium, 40°, 3 step design, green



G-24314 PEP 2.4 MM INFUSION SLEEVE
Pure Efficiency Phaco
for 2.4 mm incisions, silicone



G-34089 SINGLE-USE 2.4 MM PHACO KNIFE
Pure Efficiency Phaco
angled, 2.4 mm, bevel up
6 pcs. per box, sterile



G-34176 SINGLE-USE 19G MVR STILETTO
double bevel, 19 Gauge / 1.1 mm
6 pcs. per box, sterile



G-34181 SINGLE-USE 20G MVR STILETTO
double bevel, 20 Gauge / 0.9 mm
6 pcs. per box, sterile



G-24332 PEP 2.4 MM I/A HANDPIECE
Pure Efficiency Phaco
for 2.4 mm incisions, straight design
with G-24314 PEP infusion sleeve



G-24432 PEP 2.4 MM I/A HANDPIECE
Pure Efficiency Phaco
for 2.4 mm incisions, angled design
with G-24314 PEP infusion sleeve

for 2.2 mm + 2.4 mm incisions

G-24110 WRENCH for reusable ultrasonic tips

G-24122 TEST ANTERIOR CHAMBER silicone

G-28431 ULTRASONIC HANDPIECE P3 B diameter 16 mm, length 133 mm, large plug

G-30450 ULTRASONIC HANDPIECE ERGO for MEGATRON, MEGATRON S3 and MEGATRON S4 HPS systems,
diameter 12 mm, length 129 mm, weight 42 g, large plug

GEUDER AG reserves the right to make changes to technical details in response to recent developments. Geuder does not assume liability for the accuracy of each individual statement.
Illustrations not drawn to scale.